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Application and safety instructions for the MEISINGER

EN

# Benex®-Control

## Root Extraction System

Art.-No.: CBE00

Benex®-Control is outstandingly suitable for the controlled and safe extraction of dental roots and dental root fragments, so the conventional extraction method using forceps and leverage can be avoided. Thanks to the innovative construction of the extractor, the root can be easily removed in a particularly controlled manner without damaging bone or soft tissue. To guarantee an optimal vertical transmission of force, there must be a flat area providing a supporting surface to the extractor. The quadrant support was developed specifically to provide this optimal support while also protecting existing dentition and providing support for the extractor in case of edentulous circumstances.

## Indications


The Bone Management® system Benex®-Control is indicated for the extraction of single-root tooth roots and tooth root fragments. In the case of multi-root tooth roots, these are to be separated before extraction.

## Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.


### Preparation of the tooth root

The desmodontal fibers (Sharpey's fibers) in the sulcus are cut circularly with the aid of a periotome. Strong, long roots have to be loosened or luxated axially in both mesial and distal directions for a period of 30 seconds using a delicate lifting device / twister. In this process, care must be taken that no transversal movements are made. When extracting molars, the roots have to be separated before removal.


### Placing of an axial drill hole into the tooth root

With the aid of a fitting pilot bur (A2001 for the Extraction screws BE001230 060 and BE001230 070, B2001 for the Extraction screws BE001230 080 and BE001230 085) an axial drill hole has to be drilled into the center of the dental tooth root or tooth root fragment. This should have a depth of about 7 mm in hard tissue. A deeper drill hole is not necessary. The use is carried out with sufficient cooling, preferably with a physiological, sterile saline solution. In order to remove the drilling chips easily, the pilot bur has to be used intermittently, i.e. with up and down movements.



### Insertion of the extraction screw

Depending on the case situation, the short or long extraction screw with the appropriate diameter is selected and can be carefully inserted into the tooth root using the jollying aid. Note: Care must be taken that the jollying aid is pressed onto the extraction screw. This is the only way to ensure that there is a sufficient level of clamping force for safe transport into the oral cavity.


### Connection of the traction string with the protection cap

The images illustrate how the traction string is threaded into the protection cap. For ease of handling, it is an advantage to keep the traction string very short.


### Connection of the traction string and the protection cap with the extraction screw

The protection cap is moved slightly upwards. Then the traction string is hooked into the extraction screw. To secure the connection, the protection cap is now pushed over the screw head.


### Positioning of the extractor

The extractor is positioned on the adjacent tooth crowns with its tooth-friendly Teflon® surface. The opening of the round, rotating segment plate is oriented vestibularly, so the extraction screw should be clearly visible. The traction string is led through the pulley and attached to the extraction slide hook. After final positioning, the screw and the traction string must have the same axial direction. By turning the hand screw, the traction string is put under light pull, so that it does not hang out.


### Extraction of the tooth root

Once the extractor has been correctly positioned, extraction takes place by slow, controlled turning of the hand screw. In particular, in the case of long, strong roots it is important to preload the periodontal fibers with a submaximal pulling force at intervals of 30-40 seconds. The extractor is positioned so that it is supported under tension. Hold this position until the periodontal fibers give in after about 30-40 seconds. Subsequently, the periodontal fibers are preloaded by further rotation of the hand screw. This process is repeated until the tooth root has been completely extracted.

### Optional application of the quadrant support in case of an uneven supporting surface

In order to guarantee an optimal vertical transmission of force, there must be a right-angled surface for the application of the extractor. In case of a bad and uneven supporting surface or a free-end situation or an edentulous jaw, the quadrant support provides the optimization of the right-angled supporting surface. It allows for the even distribution of force to protect the teeth. The use of the quadrant support is often necessary for extractions in the lower jaw and palatal roots in the upper jaw.

### Utilization

- Fill the quadrant support with a fast-setting silicon or alginate and model a gap in the area of the buccal hole region.
- Position the quadrant support with the gap over the root to be extracted, right-angled to the root axis onto the jaw.
- Remove the quadrant support when the silicon or alginate starts setting. The entire setting process takes place extraorally.
- After root luxation and the following insertion of the extraction screw, the repositioning of the quadrant support and the positioning of the extractor follow. The tooth extraction is done according to the procedure in the point „Extraction of the tooth root“.

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
# Instructions for disassembly and assembly for the processing of the Benex extractor

For optimal processing, the Benex Extractor must be completely disassembled prior to cleaning, disinfection and sterilization. To do this, first remove the knurled screw from the extractor by turning it counterclockwise off the thread.


After that, the rotary segment can be removed from the main body by turning it clockwise or counterclockwise. A certain direction of rotation is not specified.

By pulling or pushing the two cross-pins are to be detached from the body of the Benex extractor. It should be noted that the guide cover of the threaded pin has to be secured during this procedure to prevent the unwanted expansion of the compression spring.


Now the guide cover can be removed carefully. The tensioning element should be fixed. After removing the compression spring, the tensioning element and the two O-rings can be removed.

Subsequently, the Teflon segment plate can be released from the rotary segment by gently pulling.

Now all individual parts of the extractor can be cleaned and disinfected according to valid processing instructions. After successful cleaning and disinfection and prior to packaging and sterilization of the Benex Extractor, the individual parts can be assembled. During assembly, it must be ensured that the O-rings have to be inserted prior to the cross pins in order to guarantee a perfect fit.

When preparing the Benex-Control Extractor for sterilization, care must be taken to ensure that the ball does not engage in a depression of the segment plate but is in a semi-oblique position.

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## Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

### Absolute contraindications

- Acute vascular diseases e.g. Myocardial infarction, severe heart failure, stroke
- Hematopoietic diseases such as acute leukemia and agranulocytosis
- Radiation and chemotherapy in the head and neck area
- Use in the case of bone fractures
- Dental ankylosis

### Relative contraindications

- Poor general condition of the patient (e.g. anemia)
- Acute infectious diseases as well as acute oral cavity infections
- Congenital tendency to for infections (e.g. hypimmunity, immunosuppressive therapy)
- coagulopathy, anticoagulant therapy
- Poor condition and degree of damage of the adjacent teeth (e.g. loosening, unstable restorations)

### Special instructions for processing

During pre-cleaning and manual cleaning and disinfection the Benex-Control traction strings must be bent several times in all directions.

### General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

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