Absolute Kontraindikationen Absolute contraindications Contraindicaciones absolutas:

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• Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)

• Active infections as well as pathological processes

• Insufficient bone supply (quality / quantity) of the donor site as well as insufficient bone quality of the donor region

• Diseases that impair bone metabolism

IMPORTANT: Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent tooth / roots (risk of damage, infection / dehiscences)

NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

CAUTION: For all cutting and drilling procedures (transplant removal / other work), the following applies: In order to reduce the risk of overheating the bone and necrosis formation along with it, the respective instrument has to be worked with intermittently using little pressure and subject to continuous cooling using a sterile physiological saline solution. Trepanning speed may not exceed 300 min⁻¹.

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
¹204=RA, 205=RA L ²Largest working part diameter in 1/10 mm Internal diameter External diameter

Lagerfräser Ablative burs Fresas ablativas					Diamantierte Sägen Diamond saws Sierras diamantadas			Bohrer Verlängerung Drill extension Extensión de fresa		Schraubendreher Screw driver Atornillador
Fig.	TC21X	TC21X	TC21X	TC21X	TC21X	231DC (*)	231DC (*)	231DC (*)	BV028	SW150
Shank ¹	205	205	205	205	205	204	204	204		
Size ²	040	050	060	070	080	070	100	130		
Length mm	14.0	14.0	14.0	14.0	14.0	0.3	0.3	0.3	28.0	

¹204=RA, 205=RA L ²Largest working part diameter in 1/10 mm

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¹ 204=RA, 205=RA L ²Largest working part diameter in 1/10 mm

Vorkörner Initial bur Fresa inicial		Trepan Trepine Trépano	Lagerfräser Ablative bur Fresa ablativa
Fig.	229VT	229	207RF
Shank ¹	205	205	205
Size ²	075	075	075
Length 1 (mm)	10.0	10.0	3.2
Length 2 (mm)	2.2		3.0
	7.5	7.5	
	8.5	8.5	7.5

¹205=RA L ²Largest working part diameter in 1/10 mm Internal diameter External diameter

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- Secure and easy transplantation of autogenous bone cylinders with a set of precise tools
- Optimally suitable for horizontal and vertical augmentative grafts
- Precise and predictable application due to depth marks
- Suitable for preparing the bone for subsequent insertion of all standard forms of intraoral dental implant fixtures
- Simple, safe transplant of bone cylinders for the restorative construction of a vertically atrophied jawbone using intelligently coordinated instruments
- Prepare a fresh flat contact surface and a central insertion hole for the osteosynthesis screw in a single operation
- Bone cylinder is securely located with no gap
- Optimum site conditions for secure healing
- Considerably reduced operation time due to systematic procedure
- Optimum combination with Bone Management® set Screw System TX

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany


Hager & Meisinger GmbH | Hansemanstr. 10 | 41468 Neuss | Germany
Tel.: +49 (0) 21 31 20 12-0 | Fax: +49 (0) 21 31 20 12-222 | www.meisinger.de | info@meisinger.de

Meisinger USA, L.L.C. | 10200 E. Easter Avenue | Centennial, Colorado 80112 | USA
Tel.: +1 (303) 268-5400 | Toll free: +1 (866) 634-7464 | Fax: +1 (303) 268-5407
www.meisingerusa.com | info@meisingerusa.com

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Transfer-Control Plus is a combination of the approved Transfer-Control and Transfer-Ring-Control systems. Instruments that are perfectly matched with each other make perfectly fitting bone transplants possible and provide the best conditions for quick and safe healing by creating a congruent and fresh recipient bone site. Perfectly transplanted bone cylinders of this type result in a bone site that is capable of implantation after only approx. 3–4 months due to quicker vitalization and healing. Transfer-Control Plus contains a variety of different instrument sizes for a high level of flexibility when extracting individually adjusted transplants.

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
The Transfer-Control system is indicated for the intraoral extraction of bone cylinders as well as for horizontal and vertical bone deficits in the maxilla and the mandible.

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

In cases of severe horizontal alveolar ridge resorption (where even bone-spreading techniques cannot be deployed) **horizontal augmentative grafts** are required for compensation purposes.

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
Depending on the implant diameter required, initially a bone cylinder of the correct diameter is extracted using an appropriate trephine, e.g. from the retromolar region of the mandible. If required, this can be loosened from the donor bone with the osteotomy disc or the osteotomy diamond bur. NOTE: Should no initial bur be used, the trephine should be initially used rotating toward the left (fewer rotations). Trepanning speed may only be 300 min⁻¹.

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
Next, using the ablative bur corresponding to the trephine used previously, a bone site is ground on the horizontally resorbed alveolar process at the tooth axis. The surgical technique and the properties of the ablative bone bur lead to a hyperaemic environment, which accelerates the wound healing in post transplantation. The donor site may be filled with collagen or other grafting material. The intervention is terminated by a current gingival suture.

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
The extracted bone cylinder is placed into the recipient bone site to compensate the horizontal bone deficit.

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
Depending on the case situation at hand, the bone cylinder can be attached with a fixation screw from the Screw System TX (Art.-No. BTX00, BBTXPR), which optimally compliments Transfer-Control.

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
In cases of severe vertical alveolar ridge resorption, **vertical augmentative grafts** are required for compensation. In this case the appropriate wheel-cutter creates a congruent cavity in the crestal alveolar ridge. One should take note of the correlating trephines to be used. Prepare an autogenous bone graft in the same manner as described for horizontal defects.

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
The Transfer-Ring-Control system is indicated for intraoral bone cylinder transplantation in the case of vertical bone defects.

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
The initial phase requires the donor region to be marked with the initial bur trephine. The pilot drilling thus created is used later to ensure perfect centering of the osteosynthesis screw in the implant site. At the same time, the surface of the intended bone cylinder is ground flat thus creating a radial groove.

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
Next, the trephine is used to cut a bone cylinder of the correct diameter from the donor site. Depending on the vertical requirements at the site, this cylinder can have a height of 6 mm to 14 mm. The diameter of the bone cylinder can range from 4 mm up to 8 mm. Trepanning speed may a maximum of 300 min⁻¹.

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With the help of the hand instrument, the bone cylinder is carefully taken from the donor region, and then subsequently the wound area is treated if necessary.

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
By using the end milling cutter, surface roughness of the bone cylinder can be milled flat

in a controlled manner accordingly. If a fixation of the bone cylinder is necessary, the next step is to drill a central guide hole for the osteosynthesis screw.

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
Next, the end milling cutter is used to form a standard sized bone site with a flat, fresh contact area into which the bone cylinder can be inserted perfectly. Since the

instruments used are all designed to complement each other, the transplant will be a perfect, gap-free fit.

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After inserting the bone cylinder into the recipient bone site, this is attached using an appropriate osteosynthesis screw from the Screw System TX (Art.-No. BTX00, BBTXPR), which optimally compliments Transfer-

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If required, splinters on the upper surface of the bone cylinder can be removed using the diamond polisher or the allport bur which are included in the set.

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